

BETHYL | **FORTIS**
LABORATORIES, INC. | LIFE SCIENCES

GeneID 280717

APPLICATIONS	IEP, DD
SPECIES REACTIVITY	Bovine
AMOUNT	2 ml
STORAGE/SHELF LIFE	2 – 8°C / 2 years from date of receipt
PHYSICAL STATE	Liquid
BUFFER	Serum containing 0.09% Sodium Azide
ISOTYPE	IgG
ORIGIN	USA
PRODUCTION PROCEDURES	By immunoelectrophoresis, the antiserum reacts specifically with bovine albumin. No antibody was detected against other bovine serum proteins. This antibody may cross react with albumin from other species.
APPLICATIONS	Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use. IEP Neat DD Neat
APPLICATION NOTES	For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion Ouchterlony.
ADDITIONAL INFO	https://www.fortislife.com/p/A10-127 Use the link above to view SDS, a current list of citations, and other product specific information.

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD Date: August 29, 2023