

Sheep IgG Heavy and Light Chain Antibody

Donkey Polyclonal

Antiserum

Catalog No. A130-100

Lot No. 241119

| | | | | | |
|------------------------------|--|-----|------|----|------|
| APPLICATIONS | IEP, DD | | | | |
| SPECIES REACTIVITY | Sheep | | | | |
| AMOUNT | 2 ml | | | | |
| STORAGE/SHELF LIFE | 2 – 8°C / 2 years from date of receipt | | | | |
| PHYSICAL STATE | Liquid | | | | |
| BUFFER | Serum containing 0.09% Sodium Azide | | | | |
| ISOTYPE | IgG | | | | |
| ORIGIN | USA | | | | |
| PRODUCTION PROCEDURES | <p>By immunoelectrophoresis, the antiserum reacts specifically with Sheep IgG and with light chains common to other rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.</p> <p>This antibody may cross react with IgG from other species.</p> | | | | |
| APPLICATIONS | <p>Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.</p> <table><tr><td>IEP</td><td>Neat</td></tr><tr><td>DD</td><td>Neat</td></tr></table> | IEP | Neat | DD | Neat |
| IEP | Neat | | | | |
| DD | Neat | | | | |
| APPLICATION NOTES | For use in precipitin gel reactions (e.g. immunoelectrophoresis (IEP) or double diffusion Ouchterlony (DD)). | | | | |
| ADDITIONAL INFO | <p>Optimal working dilutions should be determined experimentally by the investigator.</p> <p>https://www.fortislife.com/p/A130-100</p> <p>Use the link above to view SDS, a current list of citations, and other product specific information.</p> | | | | |

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.
Michael Spencer, PhD

Date: November 22, 2024