

Sheep IgG–Fab Antibody

Rabbit Polyclonal

Antiserum

Catalog No. A130-107

Lot No. A130-107-2



APPLICATIONS	IEP, DD				
SPECIES REACTIVITY	Sheep				
ISOTYPE	IgG				
AMOUNT	2 ml				
STORAGE/SHELF LIFE	2 – 8° C / 2 years from date of receipt				
PHYSICAL STATE	Liquid				
BUFFER	Serum containing 0.09% Sodium Azide				
ORIGIN	USA				
PRODUCTION PROCEDURES	<p>Antiserum was solid phase adsorbed to ensure specificity.</p> <p>By immunoelectrophoresis, the antibody reacts specifically with Fab fragment of sheep IgG and with light chains common to other sheep immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins.</p> <p>This antibody may cross react with Fab Fragments from other species.</p>				
APPLICATIONS	<p>Centrifuge tube to remove product from lid. Optimal working dilutions should be determined experimentally by the investigator. Prepare working dilution immediately before use.</p> <table><tr><td>IEP</td><td>Neat</td></tr><tr><td>DD</td><td>Neat</td></tr></table>	IEP	Neat	DD	Neat
IEP	Neat				
DD	Neat				
APPLICATION NOTES	For use in precipitin gel reactions; such as, immunoelectrophoresis and double diffusion Ouchterlony.				
ADDITIONAL INFO	<p>https://www.bethyl.com/product/A130-107</p> <p>Use the link above to view SDS, a current list of citations, and other product specific information.</p>				

This document certifies that this product has met all of the quality control standards defined by Bethyl Laboratories, Inc.

Eric McIntush, PhD | Chief Scientific Officer

Date: December 3, 2018