



Detection of human HDMX by western blot of immunoprecipitates. *Samples:* Whole cell lysate (1 mg for IP; 20% of IP loaded) from Jurkat cells. *Antibodies:* Affinity purified goat anti-HDMX antibody A303-870A (lot A303-870A-1) used for IP at 6 μ g/mg lysate. HDMX was also immunoprecipitated by rabbit anti-HDMX antibody A300-287A For blotting immunoprecipitated HDMX, A300-287A was used at 1 μ g/ml. *Detection:* Chemiluminescence with an exposure time of 3 minutes.