

AbNano[®] VHH Anti-NK-Cell Library

Catalog No. L770-200

APPLICATIONS	Phage display affinity selections for VHH Discovery
AMOUNT	1 mL
CONCENTRATION	1 x 10 ¹² plaque-forming units (PFU) per ml
STORAGE/SHELF LIFE	Frozen at -80 °C / 1 year from date of receipt
PHYSICAL STATE	Liquid
BUFFER	Phosphate-buffered saline plus glycerol
ISOTYPE	VHH from llama and alpaca IgG2/IgG3
ORIGIN	USA
PRODUCTION PROCEDURES	Variable domains were isolated recombinantly and transformed as phagemid library. Phagemid were infected with M13KO7 filamentous bacteriophage. Secreted bacteriophage cultures were PEG-precipitated to high concentration, then normalized to 1E12 pfu/mL by titer dot assay. Centrifuge tube to remove product from lid. 2.5E10 pfu or more is advised for input panning rounds.

For in vitro laboratory use only. Not for any clinical, therapeutic or diagnostic use in humans or animals. Not for human or animal consumption. This product may not be resold or modified for resale without the prior written approval of Bethyl Laboratories, Inc. The information provided in this data sheet is believed to be correct but does not purport to be all-inclusive and is intended to be used as a guide. Bethyl Laboratories, Inc. shall not be liable or responsible in any way for use of either this information or the material supplied. Disposal of hazardous material may be subject to federal, state or local laws or regulations.